



METABOLITLISTA LITA URINTESTER

Analyter	Koncentration (ng/mL)	Läkemedel som ger positivt utslag
AMFETAMIN (AMP 300)		
D,L-Amfetaminsulfat	75	Attentin
L-Amfetamin	10,000	
(±) 3,4-Metylendioxyamfetamin	150	
Fentermin	300	
Maprotilin	15,000	(Ludiomil)
Metoxyfenamin	2,000	
Dexamfetamin (Lisdexamfetamin)	300	Elvanse
BENSODIAZEPINER (BZO 200)		
Alprazolam	70	Alprazolam, Xanor
a-hydroxyalprazolam	1,000	
Klobazam	120	(Frisium; licenspreparat)
Klonazepam	300	Iktorivil-preparat enbart, (spåras ej via urinen)
Klorazepat-dipotassium	300	
Delorazepam	600	
Desalkylflurazepam	120	
Flunitrazepam	120	Flunitrazepam
(±) Lorazepam	2,000	Lorazepam, Temesta
RS-Lorazepamglukuronid	120	
Midazolam	4,000	(Midazolam, BUCCOLAM)
Bromazepam	600	
Klordiazepoxid	600	
Nitrazepam	120	Nitrazepam, Apodorm, Mogadon
Norklordiazepoxid	70	
Nordiazepam	600	Diazepam, Stesolid
Oxazepam	200	Oxascand, Sobril
Temazepam	70	
Diazepam	200	Diazepam, Stesolid
Estazolam	4,000	
Triazolam	2,000	Halcion

BUPRENORFIN (BUP 10)		
Buprenorfin	10	Buprenorfin, Norspan, Suboxone, Temgesic, Subutex, Bupremyl, Buprefarm, Bunalict, Espranor, Zubsolv
Buprenorfin 3-D-Glukuronid	50	
NorBuprenorfin	50	
NorBuprenorfin 3-D-Glukuronid	100	
KOKAIN (COC 100)		
Bensoylekgonin	100	
Kokain HCl	80	
Komaetylen	7,000	
Ekgonin	10,000	
MARIJUANA (THC 200)		
Cannabinol	140,000	
11-nor- Δ 8-THC-9 COOH	120	
11-nor- Δ 9-THC-9 COOH	200	Sativex, Triumeq
Δ 8-THC	68,000	
Δ 9-THC	68,000	
Efavirenz	+	Efavirenz, Atripla, Stocrin, Emtenef, Padviram
Cannabidiol (CBD-olja)	+	Epidyolex; kan korsreagera då en del av CBD metaboliseras till en strukturell analog av THC.
MARIJUANA (THC 50)		
Cannabinol	35,000	
11-nor- Δ 8-THC-9 COOH	30	
11-nor- Δ 9-THC-9 COOH	50	Sativex, Triumeq
Δ 8-THC	17,000	
Δ 9-THC	17,000	
Efavirenz	+	Efavirenz, Atripla, Stocrin, Emtenef, Padviram
Cannabidiol (CBD-olja)	+	Epidyolex; kan korsreagera då en del av CBD metaboliseras till en strukturell analog av THC.
MARIJUANA (THC 25)		
Cannabinol	17,500	
11-nor- Δ 8-THC-9 COOH	15	
11-nor- Δ 9-THC-9 COOH	25	Sativex, Triumeq

Δ8-THC	8,500	
Δ9-THC	8,500	
Efavirenz	+	Efavirenz, Atripla, Stocrin, Emtenef, Padviram
Cannabidiol (CBD-olja)	+	Epidyolex; kan korsreagera då en del av CBD metaboliseras till en strukturell analog av THC.

METAMFETAMIN (MET 300)		
ρ-Hydroxymetamfetamin	7,500	
D-Metamfetamin	300	
L-Metamfetamin	6,000	
(±)-3,4-Metylendioxy-Metamfetamin	3,750	
Mefentermin	15,000	
Dexamfetamin (Lisdexamfetamin)	+	Elvanse
MDMA 500 (Ecstasy)		
(±) 3,4-Metylendioxy Metamfetamin HCl	500	
(±)3,4-Metylendioxy Amfetamin HCl	3,000	
3,4-Metylendioxyetyl-Amfetamin	300	
OPIATER (OPI 200)		
Kodein	160	Kodein, Altermol, Ardinex, Citodon (Forte), Panocod, Spasmofen, Treo Comp, Paracetamol/Kodein Evolan
Levorfanol	1,000	
Morfin-3-β-D-Glukuronid	600	
EtylMorfin	4,000	Cocillana-Etyfin, Lepheton
Hydrokodon	40,000	
Hydromorfon	2,000	Palladon
6-Monoacetylmorfin	200	
Norkodein	4,000	
Normorfon	40,000	
Oxikodon	20,000	(Oxikodon, Oxycodone, OxyContin, OxyNorm, Targiniq, Armoneve)
Oxymorfon	40,000	
Prokain	10,000	
Tebain	4,000	
Morfin	200	Morfin, Oramorph, Dolcontin, Depolan, Spasmofen

TRAMADOL (TML 100)		
n-Desmetyl-cis-tramadol	200	
Cis-tramadol	100	Tradolan, Tiparol, Nobligan, Gemadol, Dolatramyl
Procyklidin	100,000	
o-Desmetyl-cis-tramadol	10,000	
Fencyklidin	100,000	
d,l-O-Desmetyl venlafaxin	50,000	Efexor, Venlafaxin
Fexofenadin	3,900	Allegra, Altifex, Fexofenadin, Nefoxef, Telfast
KETAMIN (KET 1000)		
Ketamin	1,000	Ketalar, Ketanest
Dextrometorfan	2,000	
Metoxyfenamin	25,000	
d-Norpropoxifen	25,000	
Promazin	25,000	
Prometazin	25,000	Lergigan
Pentazocin	25,000	
Fencyklidin	25,000	
Tetrahydrosolin	500	
Mefentermin	25,000	
(1R, 2S) - (-)-Efedrin (Efedrinhydroklorid)	100,000	Efedrin, Lepheton, Mollipect
Disopyramid	25,000	Durbis Retard
Bensfetamin	25,000	
(+) Klorfeniramin	25,000	
Klonidin	100,000	Catapresan
EDDP	50,000	Metadon
4-Hydroxyfencyklidin	50,000	
Levorfanol	50,000	
MDE	50,000	
Meperidin	25,000	
d-Metamfetamin	50,000	
l-Metamfetamin	50,000	
3,4- Metylendioxymetamfetamin (MDMA)	100,000	
Tioridasin	50,000	

OXIKODON (OXY 100)		
Oxikodon	100	Oxikodon, Oxycodone, OxyContin, OxyNorm, Targiniq, Armoneve
Oximorfon	300	
Levorfanol	50,000	
Hydrokodon	25,000	
Hydromorfon	50,000	Palladon
Naloxon	25,000	Naloxon, Nexodal, Suboxone, Targiniq, Zubsolv, Prenoxad, Oxycodone/Naloxone Sandoz, Nyxoid, Bunalict, Armoneve
Naltrexon	25,000	Naltrexon, Mysimba
EDDP 100		
2-Etylidin-1,5-dimetyl-3,3-difenylpyrrolidin (EDDP)	100	
Prometazin	+	Lergigan
FENTANYL (FYL 10)		
Alfentanil	300,000	Alfentanil, Rapifen
Fenfluramin	25,000	
Norfentanyl	10	
Buspiron	8,000	Buspiron, Anksilon
Fentanyl	50	Fentanyl, Abstral, Actiq, Durogesic, Effentora, Instanyl, Leptanal, Matrifen, PecFent
Sufentanil	25,000	Sufenta
Paliperidon	10,000	Invega, Trevicta, Xeplion, Paliperidon
Fexofenadin	100,000	Allegra, Altifex, Fexofenadin, Nefoxef, Telfast
Risperidon	500	Risperidon, Rispemyl, Risperdal, Rispolept
Prometazin	≥1 250 000	Lergigan, Prometazin
Ciprofloxacin	+	Ciprofloxacin, Ciproxin, Cetraxal, Ciloxan
Klozapin	+	Clozapine, Froidir, Leponex
Perfenazin	+	Trilafon
Metylfenidat	3,000	Concerta, Equasym, Medanef, Medikinet, Metylfenidat, Metynor, Ritalin
Bromhexin	75,000	Bisolvon, Bromhexin, Coldmexin, Coldypect, Etrilect, Mollipsect, Mollisupradin
Etyl- β-D-Glukuronid (ETG 500)		
Etyl- β -D-Glukuronid	500	
Morfin 3β-glukuronid	100,000	
Glucuronic Acid	100,000	
Metanol	>100,000	

Propyl β -D-glukuronid	50,000	
Morfin 6β -glukuronid	100,000	
Etanol	>100,000	

KLONAZEPAM (ACL 100)

a-hydroxyalprazolam	2000	
Bromazepam	2000	
Klordiazepoxid	2000	
Klobazam	3000	
Klonazepam	800	
Delorazepam	2000	
Desalkylflurazepam	2000	
Flunitrazepam	1000	Flunitrazepam
Lamotrigin	+	Lamotrigin
Naproxen	+	Naproxen
Norklordiazepoxid	1500	
Nordiazepam	5000	
RS-Lorazepam	900	Lorazepam, Temesta
Temazepam	3000	
7-aminoklonazepam	100	Iktorivil
7-aminoflunitrazepam (f.d. Rohypnol)	10,000	Flunitrazepam

LSD

LYSERGSYRADIETYLAMID	10	
Buspiron	+	Buspiron

METYLFENIDAT (MPD 150)

Metylfenidat	150	Metylfenidat, Concerta, Ritalina, Medikinet, Equasym Depot, Medanef
Ritalinsyra	1,000	

ZOLPIDEM (ZOL 50)

Zolpidem	50	Zolpidem, (Edluar), Stilnoct
----------	----	------------------------------

SPICE/K2

<i>Se produktblad för Spicesorterna</i>		
Lamotrigin	Okänd nivå	Lamotrigin, Lamictal, Keppra

ZOPIKLON (ZOP 50)

Zopiklon	50	Imovane, Zopiklon/Zopiclone
----------	----	-----------------------------

PREGABALIN (PGB 500ng/mL)

Pregabalin	500	Lyrica, Pregabalin, Brieka
Gabapentin	20 000	Gabapentin, Neurontin

GABAPENTIN 2000 ng/mL

Gabapentin	2000	Gabapentin, Neurontin
Pregabalin	<100 000	Lyrica, Pregabalin, Brieka
Fenfluramin	>10 000	
Venlafaxin	>100 000	Efastad, Efexor, Venlafaxin
Fenetylamin	>100 000	
Tetrakain	>500 000	Rapydan, Mekostest, Tetrakain
Amoxicillin	>10 000 000	Amimox, Amoxicillin, Betaklav, Bioclavid, Imacillin, Klaximol, Nexium HP, Spektramox
Isoxsuprin	1 000 000	

KAT 150 ng/mL

Norpseudoefedrin Hcl	150	
p-Hydroxyamfetamin	100	
3,4-Metylendixoyamfetamin	100	
Tryptamin	12 500	
Metoxyfenamin	12 500	
D.L-Amfetamin	100	Elvanse, Attentin, Metamina

GHB (OBS! Kemiskt test med färgtolkning)

Alkohol		
Vitamin C		

PCP

4-Hydroxyfencyklidin	12,500	
Fencyklidin	25	